STERILE GOODS  
MANAGEMENT

# AESCULAP<sup>®</sup> BipoJet<sup>®</sup>

BIPOLAR INSTRUMENTS

AESFULAP® BIPOLAR INSTRUMENTS  
AESFULAP® BIPOLARE INSTRUMENTE

BipoJet® bipolar instruments represent technical solutions and improved product characteristics, compared to standard bipolar instruments, for the benefit of the user.

Particularly important are the ceramic components used in all BipoJet® instruments. Ceramic is a hard and heat resistant material.

For this reason, we use special ceramic components for our BipoJet® instruments to insulate the working area. For forceps, the practical advantages are the lower thermal spread in the coagulation zone and the reduction of collateral coagulation damage compared to conventional bipolar instruments<sup>1</sup>. BipoJet® scissors have an insulation to allow simultaneous dissection and coagulation.

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

#### Defined coagulation area – lower thermal spread<sup>1</sup>

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Conventional Forceps  
Herkömmliche Pinzette

BipoJet® Forceps  
BipoJet® Pinzette

<sup>1</sup> Both figures were generated in vitro on animal tissue with the the generator GN640 setting 10W / 15W, activation 8 sec. The ceramic insulated forceps showed a lower temperature peak as the uninsulated forceps.  
Insulated: 10W: ~80°C, 15W: ~95°C  
Uninsulated: 10W: ~98°C, 15W: ~180°C

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## AESCULAP® BIPOLARE INSTRUMENTE

### BIPOJET® DISSECTING SCISSORS | BIPOJET® PRÄPARIERSCHEREN

#### Cutting and coagulation in one step

#### Without BipoJet®:

Example of daily OR routine with three instruments (dissecting scissors, dissecting forceps\*, monopolar HF electrode).

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#### With BipoJet®:

Just one instrument for dissecting and coagulation.

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#### Anti-twist cable

The BipoJet® cable and jack are designed and constructed in such a way that the cables cannot become twisted. This means that the scissors can still be opened even after being rotated repeatedly.

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
\* Not the intended use of non-insulated dissecting forceps – risk of collateral coagulation.

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## Parallel cutting and coagulating

Parallel cutting and coagulation of the incised tissue: The tissue is cut while coagulation mode of the HF generator is activated.

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
## Precoagulation

Prior to performance of the cut, the scissors are rotated so as to coagulate e. g. more robust tissue and larger blood vessels.

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
## Pinpoint coagulation

Bleeding can be precisely controlled in localized areas with the scissors tip closed or slightly open.

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
## Surface coagulation

Surface coagulation performed between the open blades of the scissors.

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#### **BC691R**

##### **METZENBAUM**

extra fine tip, 120 mm, 4 3/4"  
extra feine Spitze, 120 mm

#### **BC692R**

##### **METZENBAUM**

extra fine tip, 140 mm, 5 1/2"  
extra feine Spitze, 140 mm

#### **BC693R**

##### **METZENBAUM**

extra fine tip, 160 mm, 6 1/4"  
extra feine Spitze, 160 mm

#### **BC694R**

##### **STEVENS**

Tenotomy tip, 120 mm, 4 3/4"  
Tenotomy Spitze, 120 mm

#### **BC695R**

##### **STEVENS**

Tenotomy tip, 140 mm, 5 1/2"  
Tenotomy Spitze, 140 mm

#### **BC696R**

##### **STEVENS**

Tenotomy tip, 160 mm, 6 1/4"  
Tenotomy Spitze, 160mm

**BC697R**  
METZENBAUM  
extra fine tip, 180 mm, 7"  
extra feine Spitze, 180 mm

**BC698R**  
METZENBAUM  
extra fine tip, 230 mm, 9"  
extra feine Spitze, 230 mm

**BC699R**  
METZENBAUM  
extra fine tip, 280 mm, 11"  
extra feine Spitze, 280 mm

# AESCULAP® BIPOLAR INSTRUMENTS

## AESCULAP® BIPOLARE INSTRUMENTE

### BIPOJET® DISSECTING SCISSORS | BIPOJET® PRÄPARIERSCHEREN

#### **BC685R**

METZENBAUM  
standard tip, 180 mm, 7"  
standard Spitze, 180 mm

#### **BC687R**

METZENBAUM  
standard tip, 230 mm, 9"  
standard Spitze, 230 mm

#### **BC689R**

METZENBAUM  
standard tip, 280 mm, 11"  
standard Spitze, 280 mm

#### **BC686R**

METZENBAUM  
fine tip, 180 mm, 7"  
feine Spitze, 180 mm

#### **BC688R**

METZENBAUM  
fine tip, 230 mm, 9"  
feine Spitze, 230 mm

#### **BC690R**

METZENBAUM  
fine tip, 280 mm, 11"  
feine Spitze, 280 mm

# AESCULAP® BIPOLAR INSTRUMENTS

## AESCULAP® BIPOLARE INSTRUMENTE

### BIPOJET® FORCEPS | BIPOJET® PINZETTEN


- Defined coagulation area.  
→ Reduction of collateral coagulation<sup>1</sup>
- Definierte Koagulationsfläche.  
→ Reduktion von kollateralen Koagulationsschäden<sup>1</sup>
- Low temperature range during coagulation.  
→ Coagulation instead of carbonization<sup>2</sup>
- Niedriger Temperaturbereich bei der Koagulation.  
→ Koagulation statt Karbonisation<sup>2</sup>

<sup>1</sup> Both figures were generated in vitro on animal tissue with the the generator GN640 setting 10W / 15W, activation 8 sec.  
The ceramic insulated forceps showed a lower temperature peak as the uninsulated forceps.  
Insulated: 10W: ~80°C, 15W~ 95°C  
Uninsulated: 10W: ~98°C, 15W~ 180°C

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
**BipoJet® vs. conventional forceps. Defined coagulation area – lower thermal spread<sup>1</sup>**

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
**Forceps with fully insulated tip.**  
Result: Confined coagulation zone.

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**Forceps with partially insulated tip.**  
Precise even in cases where the tips of the forceps have to be used for coagulation.

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
**Conventional forceps.**  
Extended coagulation zone. Undefined thermal spread.

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
**Coagulation comparison: Thermal camera images<sup>1</sup>**

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
**Ceramic forceps.** With its ceramic insulation, the BipoJet® forceps coagulates precisely at lower temperature.

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**Conventional bipolar forceps.** Clearly visible are the comparatively high temperatures at the centre, which provoke tissue adhesion, and the coagulation zone, which extends far beyond the forceps.

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\*All BiposJet® instruments are delivered with a fitting ceramic protector as standard.

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

#### AESCULAP®. QUALITY FIRST

The awareness that precision, functionality and handling of the surgical instruments produced often influence a life-saving surgical procedure is present in every one of our employees. So it is not surprising that very special instruments are created from craftsmanship and attention to detail.

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The main product trademark 'Aesculap' and the product trademark 'BipoJet®' are registered trademarks of Aesculap AG.  
Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form partial or otherwise is not permitted.

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AESCULAP® AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany