Novosyn[®] CHD

BG	Инструкции за използване	– 2	LT	Naudojimo instrukcijos	– 26
CS	Návod k použití	– 4	LV	Lietošanas instrukcija	– 28
DA	Brugsvejledning	– 6	NL	Gebruiksaanwijzing	– 30
DE	Gebrauchsanweisung	– 8	NO	Bruksinstruksjoner	– 32
EL	Οδηγίες χρήσης	– 10	PL	Instrukcja stosowania	– 34
EN	Instructions for use	– 12	PT	Instruções de uso	– 36
ES	Instrucciones de uso	– 14	RO	Instrucțiuni de utilizare	– 38
ET	Kasutusjuhend	– 16	SK	Návod na použitie	– 40
FI	Käyttöohjeet	– 18	SL	Navodila za uporabo	– 42
FR	Mode d'emploi	– 20	SV	Användarinstruktioner	– 44
HU	Használati utasítás	– 22	TR	Kullanım talimatı	– 46
IT	Istruzioni per l'uso	– 24			

CE 0123

eIFU v.01 506881 - IFU(506881 01-07/22)

B|BRAUN

SHARING EXPERTISE

B. Braun Surgical, S.A.
 Ctra de Terrassa, 121
 08191 Rubí, Spain
 Phone +34 93 586 62 00
 Fax +34 93 588 10 96
 www.bbraun.com

Novosyn® CHD

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Novosyn® CHD

Description of the device

Novosyn® CHD is a sterile multifilament braided synthetic, absorbable surgical suture material produced from a copolymer composed of 90% glycolide and 10% L-lactide (PGLA 90/10). Novosyn® CHD is coloured violet with the colouring D&C Violet No. 2 (C. I. 60725) to make it recognizable, but it is also available undyed in the natural beige colour. The braided threads are treated with an absorbable synthetic coating consisting of a mixture of equal parts of a copolymer (comprised of glycolide and L-lactide) and calcium stearate so that the suture slides easily without causing a sawing effect. Novosyn® CHD contains an antimicrobial coating of chlorhexidine diacetate at no more than 60 µg/m. Novosyn® CHD fulfils the requirements of the European Pharm. and United States Pharm. - current edition - for sterile, synthetic, absorbable sutures except for an oversize of the suture diameter.

The following table reports the average maximum oversize in diameter (mm) for Novosyn® CHD.

USP Size	EP size	Average diameter value EP (mm)	Maximum average oversize of Novosyn® CHD (mm)
5/0	1	0,100 - 0,149	0,036
4/0	1,5	0,150 - 0,199	0,041
3/0	2	0,200 - 0,249	0,036
2/0	3	0,300 - 0,349	0,031
0	3,5	0,350 - 0,399	0,056
1	4	0,400 - 0,499	0,046
2	5	0,500 - 0,599	0,031

The intended purpose of Novosyn® CHD is for soft tissue approximation and/or ligation. Novosyn® CHD is coated with the antibacterial agent chlorhexidine diacetate to reduce the risk of bacterial attachment and colonization of the suture.

The intended clinical benefits of Novosyn CHD are the following:

- Novosyn CHD provides mid-term mechanical support for at least 4 weeks, to support the physiological wound healing process.
- Novosyn CHD is coated with antibacterial agent chlorhexidine diacetate to reduce the risk of bacterial attachment and colonization of the suture. This in general would support a reduction of surgical site infections.

Indications for use

Novosyn CHD does not have a particular indication. The indications for use are related to the intended purpose of soft-tissue approximation and ligation of anatomical structures with synthetic, absorbable, braided suture materials. The patient population to be used Novosyn® CHD is adult and paediatric patients (> 2 months old) who do not meet any contraindication.

Intended user are healthcare professionals like surgeons that are qualified and familiar with the surgical techniques.

Mode of action

Suture materials are used primarily for adaptation of the wound edges to render possible an undisturbed wound healing. During the use of Novosyn® CHD sutures, a mild inflammatory reaction may occur, which is typical for an endogenous reaction to a foreign body. As time passes, the suture material is encapsulated by fibrous connective tissue. Novosyn® CHD is metabolized to glycolic acid and lactic acid by hydrolysis without causing any enduring change in the region of the wound. Approximately 75% of the original tensile strength remains after 14 days of implantation, approximately 50% after 21 days and approximately 25% after 28 days. Absorption of Novosyn® CHD is essentially complete between 56-70 days, when the tissue is normally perfused.

The effectiveness of chlorhexidine diacetate coating has been tested under *in vitro* studies and chlorhexidine diacetate in Novosyn® CHD sutures has been shown to inhibit colonization of the suture by *Staphylococcus aureus*, *Staphylococcus epidermidis*, *Methicillin Resistant S. aureus* and *Methicillin Resistant S. epidermidis* which are microorganisms known to contribute to surgical site infections. Microbiological studies carried out with guinea pigs demonstrated that surgical sutures with chlorhexidine inhibits bacterial colonization of suture after direct *in vivo* challenge with bacteria (*S. aureus*).

Contraindications

Novosyn® CHD suture materials are contraindicated for applications where prolonged support of the wound closure by the suture material is required (e.g.

cardio-vascular surgery). It is also contraindicated for ophthalmic surgery, otosurgery and also in the wound closure of cardiovascular and neurological tissues.

Novosyn® CHD suture should not be used in patients with known allergic reactions to Chlorhexidine diacetate as it could cause anaphylactic reaction. There is scientific literature advising against the use of chlorhexidine in infants aged less than 2 months, therefore Novosyn® CHD should not be used in infants aged less than 2 months.

Mode of application

To be used in accordance with surgical requirements.

Users should be familiar with the surgical procedures and techniques involving absorbable sutures when using Novosyn® CHD, as the risk of wound dehiscence may vary depending upon the site of application and the type of material used.

Adequate knot security requires the standard surgical technique of flat, square ties, with additional throws as indicated by surgical circumstances and the experience of the surgeon.

Care should be taken to avoid damage when handling surgical needles.

Grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the attachment end to the point. Grasping in the point area could impair the penetration performance and cause fracture of the needle.

Grasping at the butt or attachment end could cause bending or breakage.

Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle stick injury. In case that a take-off needle is used, after the operation is completed, please grasp securely the end of the thread and hold the needle holder, pulling strand taut, then the needle is released with a straight tug of the needle holder.

Warnings

Novosyn® CHD is a single use device. Opened unused or damaged packs should be discarded. Do not use after expiry date. Do not reuse the product: Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product.

Physicians should carefully consider the potential risks and benefits for each specific patient before using Novosyn® CHD. There are no adequate data to support the use of Novosyn® CHD in pregnant or lactating women.

No clinical studies have been performed with chlorhexidine coated sutures and preclinical studies are insufficient to claim the safety in pregnancy and during the lactation period.

Chlorhexidine diacetate is known to be ototoxic, causing damage to the neuroepithelial cells of the inner ear. Care must be taken to avoid contact of Novosyn® CHD with tissues of the inner ear.

If the needle detaches and falls into the patient and cannot be recovered, exist risk of embolism, foreign body reaction, encapsulation and potential need of additional test. Therefore, the risk is considered critical for the patient. Take-off needles should not be used for endoscopic suturing.

Precautions

Skin sutures which remain in place longer than 7 days may cause localized irritation and should be snipped off or removed as indicated.

Consideration should be taken in the use of absorbable sutures in tissues with poor blood supply as suture extrusion and delayed absorption may occur. Subcuticular sutures should be placed as deeply as possible to minimize the erythema and indurations normally associated with the absorption process.

Usage of Novosyn® CHD may not be advised in case of elderly or malnourished or debilitated patients, or in patients suffering from diseases or conditions which delay the wound healing process. When working with Novosyn® CHD suture materials, great care should be taken to ensure that the use of surgical instruments, such as forceps and needle holders, does not cause any crushing or crimping damage to the suture material.

Open the sterile pouch with caution. Avoid contact of the inner sterile part with the outer non-sterile part of the packaging.

Novosyn® CHD surgical sutures, with needle removed, contain implantable components which are considered MR safe. The implantable suture (composed of PGLA 90/10 with a coating of chlorhexidine diacetate and a coating of PGLA30/70 with calcium stearate) are non-metallic and non-conducting materials and therefore, considered MR safe.

Side Effects

As for every other suture materials, prolonged contact with salt solutions, such as urine and bile, can lead to lithiasis.

As for any other suture after implantation the following side effects may occur occasionally: transient inflammatory foreign body reaction, transient local irritation, granulation, stitch sinus, infection at the wound site or impaired

aesthetic outcomes. Existing infections may occasionally be enhanced by any foreign body. May not be excluded an occasional pain, granuloma, seroma, fibrotic scar, hardening of tissues (subcuticular tissues), irritation, hematoma, suture rejection, enhanced bacterial infectivity, wound dehiscence and haemorrhage or leakage.

In rare cases, allergic reaction, contact dermatitis or anaphylactic shock have been reported with chlorhexidine contact.

In case of contamination and/or impaired functionality of the product exists risk of injury, illness, sepsis, death and/or potential need of medical treatment or reintervention.

The following side effects may be associated with needle detachment: accidental puncture, embolism leading to a life-threatening injury, foreign body reaction, encapsulation, potential need of additional test (X-Ray) or interference with MRI scan.

Chlorhexidine diacetate is known to be ototoxic, causing damage to the neuroepithelial cells of the inner ear.

Sterilization

Novosyn® CHD is supplied sterile. Novosyn® CHD is sterilized by Ethylene Oxide and must not be resterilized. Only use Novosyn® CHD if the package is undamaged. The absence of channels along the sealing of the sterile barrier system must be checked as well as the absence of perforations on it. In case of identifying such defects, the device should be discarded in the recommended manner.

Storage

Novosyn CHD does not require any special storage and handling conditions.

Device Disposal

Once the surgical procedure is completed, the different components of Novosyn® CHD have to be discarded in special intended containers. It is the responsibility of the user to determine if disposal material is hazardous according to federal, state and local regulations. Dispose of contents and container to comply with applicable local, state, national and international regulation.

Place all needles in a sharps disposal container immediately after they have been used. Dispose of used sharps disposal containers according to your community guidelines.

Information to the User/Patient

In case that any serious incident occurs in relation to Novosyn® CHD it should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

The lifetime of Novosyn® CHD suture occurs after 56-70 days, when the complete mass of the product is absorbed.

About the qualitative & quantitative information on materials & substances which patients can be exposed are as follows: Novosyn® CHD is produced from a copolymer composed of 90% glycolide and 10% L-lactide (PGLA 90/10) and is coloured violet with the colouring D&C Violet No. 2 (0.05 – 0.20 % w/w). Novosyn® CHD is coated with a copolymer composed of 30% glycolide and 70% L-lactide (PGLA 30/70) and calcium stearate (1.0 – 5.5 % w/w) and it is also coated with chlorhexidine diacetate (at no more than 60 µg/m). Finally, the needle composition of Novosyn® CHD sutures is stainless steel AISI series 300.

Area reserved to link to the summary of safety and clinical performance on EUDAMED (Art32 MDR) [<https://ec.europa.eu/tools/EUDAMED>].

General Information

Date of information: June 2022

Symbols to be used on the label

Do not re-use

Use-by date

Batch code

Do not resterilize

Violet

Undyed

Polyglactin 90/10 and chlorhexidine diacetate

Contains a medicinal substance

Dyed/braided/coated/Absorbable

Undyed/braided/coated/Absorbable

Presence of Cobalt

Sterilized using ethylene oxide

Single sterile barrier system with protective packaging inside

Do not use if package is damaged and consult instructions for use

Date of manufacture

Consult instructions for use or consult electronic instructions for use

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Tıbbi cihaz

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Title: eIFU Novosyn CHD Initiator: Susana ? Calvillo

This document is signed electronically in compliance with the B. Braun electronic signature policies and procedures by following persons:

UserName: Calvillo, Susana (calvsues)
Title: SUT-PR-ES04-Técnico Sistemas Impresión
Date: Wednesday, 24 August 2022, 15:46 W. Europe Daylight Time
Meaning: Document signed as Author

UserName: García, Guillermo (garcgues)
Title: SU-PR-ES04-Technical support manager
Date: Tuesday, 30 August 2022, 09:08 W. Europe Daylight Time
Meaning: Precheck of Document

UserName: Orús, Silvia Salut (orussies)
Title: Global Regulatory Affairs Manager CoE CT
Date: Friday, 16 September 2022, 12:56 W. Europe Daylight Time
Meaning: Approve Document