

**B. Braun Products for manual and machine reprocessing of instruments**

Reprocessing Recommendation of the Instrument Manufacturer	Manual Reprocessing								Machine Reprocessing						
	B. Braun Product	Helipur® H plus N	Prontocid® N	Helizyme®	Cleaner N	Helipur®	Stabimed®	Helix®	Tiutol® KF	Helimatic® Neutralizer C	Helimatic® Neutralizer forte	Helimatic® Cleaner neutral	Helimatic® Cleaner enzymatic	Helimatic® Rinse neutral	Helimatic® Cleaner alkaline
acid and/or pH-neutral and/or alkaline with surfactants	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓
acid and/or pH-neutral and/or mild alkaline with surfactants	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓	✓
acid and/or pH-neutral with surfactants	✓	✓	✓	✓						✓	✓	✓	✓	✓	
pH-neutral and/or alkaline with surfactants			✓	✓	✓	✓	✓	✓				✓	✓	✓	✓
pH-neutral and/or mild alkaline with surfactants			✓	✓	✓	✓	✓	✓				✓	✓	✓	✓
pH-neutral with surfactants			✓	✓								✓	✓	✓	
Acid and/or pH-neutral and/or alkaline without surfactants								✓		✓	✓				
acid and/or pH-neutral and/or mild alkaline without surfactants										✓	✓				
acid and/or pH-neutral without surfactants										✓	✓				
pH-neutral and/or alkaline without surfactants								✓							
pH-neutral and/or mild alkaline without surfactants															
pH-neutral without surfactants															

B. Braun Medical AG  
Disinfection & Hygiene  
CH-6203 Sempach-Station

i.V. Andreas Arndt  
Head of Development and Quality Assurance

i.V. Dr. Thomas-Jörg Hennig  
Application Technology Manager Europe