



ACTREEN® HI-LITE RANGE

DISCREET & SIMPLY DESIGNED

Actreen[®] Hi-Lite Range

The invisible solution ready to stay & ready to go!

Common benefits

- Biocompatible catheter material, PVC free, DEHP free
- Non traumatic eyes for smooth bladder drainage¹
- Pre-lubricated with hydrophilic lubricant: ready to use catheter, no need for water
- Openings designed for good and reduced dexterity
- Safe handling with 2 «No touch» options.

Added value



Lightweight



Can be folded and stored inside the provided discreet pouch present inside in each box*



Convenient



2 options for opening adapted to the level of manual dexterity



No touch



2 options to avoid touching the catheter during use

Actreen® Hi-Lite Cath



To stay.

Cath



A new look

- White packaging film for more discretion



Designed for safety use

- The product can be used in aseptic conditions by removing the catheter from its packaging by the connector with sterile gloves.

Actreen® Hi-Lite Set



To go

Set



All in one concept

- Catheter + bag sealed together and folded into a compact size*



Anti reflux valve

- Integrated into the bag to avoid leakages during storage²



Graduated bag of 1 000 mL

- Approximate measurement of urine
- Easy to read graduated scale
- Bag with a freindly design, could be emptied thanks to its tear off corner⁵


*The CH18 version can not be folded.

Actreen[®] Hi-Lite Technology

A history of innovation and performance

1. A unique technology based on **hydrophobic lubricant**

Components: Water and glycerine based formula maintained throughout catheterization.



Patient's benefits

- **No water management:** gel lubricant properties
- **Ready-to-use product:** no additional water required for preparation
- **No risk of drying out before insertion³ :** the level of lubrication stays the same during the insertion

2. The catheter features and innovative packaging



Key figures


20 years
of experience


60 min
without drying³


1,3 TO 2 TIMES LESS
waste discarded/year* when using
Actreen[®] Hi Lite Products


• Multifunctional component
• Lightest product on the market*
• Design simplicity
**MAXIMUM
ECO-FRIENDLY USE**


* compared to the 3 main products used on the market


Actreen® Hi-Lite Range


Ordering information




Actreen® Hi-Lite Cath


	CH	ø (mm)	Units/box	Ref.
TIEMANN - 41 cm				
	■ 08	2,7	30 + 1 discreet pouch	238108*
	■ 10	3,3	30 + 1 discreet pouch	238110*
	□ 12	4,0	30 + 1 discreet pouch	238112*
	■ 14	4,7	30 + 1 discreet pouch	238114*
	■ 16	5,3	30 + 1 discreet pouch	238116*
	■ 18	6,0	30	238118

	CH	ø (mm)	Units/box	Ref.
NELATON - 41 cm				
	■ 08	2,7	30 + 1 discreet pouch	238208*
	■ 10	3,3	30 + 1 discreet pouch	238210*
	□ 12	4,0	30 + 1 discreet pouch	238212*
	■ 14	4,7	30 + 1 discreet pouch	238214*
	■ 16	5,3	30 + 1 discreet pouch	238216*
	■ 18	6,0	30	238218

	CH	ø (mm)	Units/box	Ref.
NELATON - 20 cm				
	■ 06	2,0	30 + 1 discreet pouch	238306*
	■ 08	2,7	30 + 1 discreet pouch	238308*
	■ 10	3,3	30 + 1 discreet pouch	238310*
	□ 12	4,0	30 + 1 discreet pouch	238312*
	■ 14	4,7	30 + 1 discreet pouch	238314*
	■ 16	5,3	30 + 1 discreet pouch	238316*

Actreen® Hi-Lite Set

	CH	ø (mm)	Units/box	Ref.
TIEMANN - 37 cm				
	■ 10	3,3	30 + 1 discreet pouch	242110*
	□ 12	4,0	30 + 1 discreet pouch	242112*
	■ 14	4,7	30 + 1 discreet pouch	242114*
	■ 16	5,3	30 + 1 discreet pouch	242116*
	■ 18	6,0	30	242118*

	CH	ø (mm)	Units/box	Ref.
NELATON - 37 cm				
	■ 10	3,3	30 + 1 discreet pouch	242210*
	□ 12	4,0	30 + 1 discreet pouch	242212*
	■ 14	4,7	30 + 1 discreet pouch	242214*
	■ 16	5,3	30 + 1 discreet pouch	242216*
	■ 18	6,0	30	242218*

* For each country, a specific letter is added at the end of the reference radical.

¹ Report n° 7.5.AA.100.1

² Report n° 7.5.AB.054.1

³ Report N° 7.5AC059.2

For safety and hygienic reasons, we advise you to use the product in the shortest possible time after opening.

This medical device is a regulated health care product and on the basis of this regulation is certified CE.

The products commercialized by B. BRAUN MEDICAL shall exclusively be used in accordance with the package inserts which have been updated and included in the boxes.

Non binding documents and photographs

Commercialized and manufactured by:

B. Braun Medical | 204 avenue du Maréchal Juin | 92100 Boulogne Billancourt - France

Tél. 01 41 10 53 00 | Fax 01 41 10 53 99 | www.bbraun.fr

B. Braun Medical SAS with a share capital of 31 000 000 € | RCS Nanterre 562050856

www.bbraun.com/continence-care-and-urology

This document, its contents, including institutional data, information, trademarks and logos mentioned herein are the exclusive property of B. Braun.

Any representation and / or reproduction, total or partial, of this document and its contents without the express prior consent of B. Braun,

is strictly prohibited and constitutes an infringement of the intellectual property rights of B. Braun.

⁴ «Results reports of the usability test product on Actreen® Hi-Lite Cath»

⁵ «Results report of the USABILITY TEST Product: Actreen® Hi Lite Set»