



**SURFACE
DISINFECTION**



Meliseptol® pure

FAST ACTING, READY TO USE SURFACE DISINFECTION

Meliseptol® pure

FAST ACTING, READY TO USE SURFACE DISINFECTION

PROPERTIES

- Highly efficient – active within 1 minute
- Enhanced cleaning performance, even on shining metallic surfaces
- Without additional antimicrobial ingredient (as aldehydes, QAC)
- Perfume-free



INDICATION

Meliseptol® pure is a highly efficient disinfectant based on a midlevel alcohol concentration for non-invasive medical inventory, medical devices e.g. examination chairs, hospital beds, theatre equipment and other kinds of surfaces. It is especially suitable for shining metallic and inaccessible surfaces. Meliseptol® pure is aldehyde- and QAC-free and contains only ethanol as antimicrobial active ingredient.

USAGE

Moisten surface completely with undiluted solution for 1 minute.

Physico-chemical data ready-to-use solution

Appearance	clear, colorless liquid
pH-value	approx. 7
Density (20 °C, g/cm ³)	ca. 0.93 g/ml
Flash point (DIN 51 755)	+ 23 °C

Meliseptol® pure – composition:

100 g solution contains: Ethanol 44.0 g, Excipients < 5 % surfactants, purified water, Ingredients in accordance with the Regulations on Detergents EC 648/2004: < 5 % non-ionic surfactants, Labelling of dangerous goods: see material safety data sheet (MSDS), Cautions: Use disinfectants safely. Always read the label and the product information before use. Do not use the product after the expiry date.

Spectrum of activity	Test Norms	Exposure Time
Bactericidal / Levurocical	EN 13697 EN 13727 EN 13624	1 min.
Mycobactericidal Tuberculocidal	EN 14348	1 min. 1 min.
Viruses (enveloped) (incl. HBV, HCV, HIV, vaccinia virus, BVDV)	DVV/RKI*	30 sec.
Viruses (non enveloped)		
▪ Rotavirus	EN 14476	30 sec.
▪ Adenovirus	EN 14476	30 sec.
▪ Norovirus (MNV)	EN 14476	30 sec.

Expert reports are available upon request

*acc. recommendation of RKI, Bundesgesundheitsbl. 01-2004, DVV/RKI suspension test

COMPATIBLE WITH THE B. BRAUN WIPES DISPENSER SYSTEM
Meliseptol® pure can be used as impregnating liquid with the B. Braun Wipes Dispenser System. The solution can remain in the dispenser system for up to 4 weeks.



Perfume-free

PRODUCT SIZE	REF
250 ml round bottle Meliseptol® pure	19926
1000 ml bottle Meliseptol® pure	19927
5 l canister Meliseptol pure	19958

B. Braun Medical AG | Infection Control | Seesatz 17 | CH-6204 Sempach

Phone +41 58 258 50 00 | Fax +41 58 258 60 00 | info.bbmc@bbraun.com | bbraun.com/infection-prevention