

# Vasco<sup>®</sup> OP eco

## Perfect fit and high quality at a value price

Anatomically shaped surgical gloves (natural rubber latex)



CoE Practitioner Use

**NEW**

**B | BRAUN**  
SHARING EXPERTISE

# Vasco® OP eco

## Perfect fit and high quality at a value price

Soft and highly flexible ergonomic latex gloves

Perfect choice for sensitivity and safety | AQL 1.0

Low chemical residual content,  
free of thiurames and mercaptobenzothiazole (MBT)



- Powder-free
- Easy donning | inner surface with latex-free polymer coating
- Safe grip | micro rough surface
- Reliable protection | 0.19 mm average wall thickness at palm
- Size-adapted glove length from 260 to 280 mm
- Sterile, packed as right/left pairs in a paper wallet and polyethylene peel pouch



size	sales unit	REF
5.5	40 pairs	6081308
6	40 pairs	6081316
6.5	40 pairs	6081324
7	40 pairs	6081332
7.5	40 pairs	6081340
8	40 pairs	6081359
8.5	40 pairs	6081367
9	40 pairs	6081375