



Vasco[®] OP Grip

FEEL THE GRIP IN DEMANDING SURGERY

Vasco® OP Grip

FEEL THE GRIP IN DEMANDING SURGERY



EXCELLENT DRY AND WET GRIP PERFORMANCE

PRECISE AND CONTROLLED INSTRUMENT GUIDANCE

LESS HAND FATIGUE DUE TO SECURE GRIP

SKIN FRIENDLY DUE TO MULTIPLE WASHING

SPECIAL INNER LINING ALLOWS EASY DONNING

YOUR CHOICE FOR SENSITIVITY AND SAFETY

- Surgical gloves according MDD 93/42/EEC, EN 455
- Protective gloves according PPE 89/686/EEC, EN 420, EN 374
- Made from natural rubber latex
- Powder-free with latex-free inner polymer coating
- Average wall thickness at palm: 0.22 mm

SIZE	SALES UNIT	REF
5.5	40 pairs	6081409
6	40 pairs	6081417
6.5	40 pairs	6081425
7	40 pairs	6081433
7.5	40 pairs	6081441
8	40 pairs	6081450
8.5	40 pairs	6081468
9	40 pairs	6081476

All Vasco® OP gloves are sterile, packed as right/left pairs in a paper wallet and polyethylene peel pouch.

CE 0123



STERILE R

CE 0321



AQL 0.65



B. Braun Melsungen AG | Carl-Braun-Str. 1 | 34212 Melsungen | Germany | www.bbraun.com

Further information on technical data, chemical and cytostatics resistance: www.bbraun.com