







# Manuplast® Vet | Technical Data Sheet

Veterinary Gloves (for single use | non-sterile)

|                              | Manuplast® Vet  | Manuplast® Vet shoulder   |
|------------------------------|---|---|
| REF                          | 6075100   | 6075151   |
| Pcs/pack                     | 100   | 50  |
|                              |   |   |
| <b>Packaging dimensions:</b> |   |   |
| Dispenser pack               | 24 x 12.5 x 6 cm (L x W x H)  | 26 x 12.6 x 6 cm (L x W x H)  |
| Transportation carton        | 36 x 25.5 x 25 cm (L x W x H)   | 32 x 26 x 28 cm (L x W x H)   |
| Transportation carton        | 10 dispenser packs  | 10 dispenser packs  |
| <b>Storage:</b>              |   |   |
| Shelf life                   | 3 years   |   |
| Conditions                   | store at room temperature, protect from dust, humidity and sun light                |   |
| <b>Glove dimensions:</b>     |   |   |
| Size                         | universal size  | universal size  |
| Length                       | 90 cm ± 1 cm  | 118 cm ± 1 cm   |
| Wall thickness               | 0.022 mm ± 0.003 mm   | 0.025 mm ± 0.003 mm   |
| <b>Glove design:</b>         |   |   |
| Shape                        | flat, ambidextrous, welded seam<br>arm-length                                       | shoulder length with neck-strap   |
|                              |  |  |
| Colour                       | transparent orange  | transparent green   |
| Surface                      | smooth  |   |



|                           |   |
|---------------------------|---|
| Glove material            | polyethylene (PE), ethylene vinyl acetate (EVA)<br><b>powder-free</b><br><b>free of latex proteins</b><br><b>free of plasticisers</b> |
| <b>Applied Standards:</b> |   |
| Tensile strength          | ISO 527-1   |
| Elongation at break       | ISO 527-3   |
| Impact resistance         | ISO 7765-1  |
| Tear resistance           | ISO 6383-2  |
| Quality certificate       | ISO 9001  |