



ProSet Safety

User Safety for CVC and Accessory Kits

ProSet Safety

Safe Instruments for Enhanced User Safety



Healthcare professionals face a higher risk of exposure to blood and body fluids in their day-to-day work. Such an exposure can lead to the transmission of pathogens, potentially resulting in an infection and causing serious damage to their health.

On May 11, 2013, the EU-Directive 2010/32/EU came into effect, which has the goal of preventing sharps injuries at hospitals and in the health sector. In Germany, the TRBA 250 guideline already requires employers to only use safe instruments in their facilities wherever it is technically possible since the year 2007. The amended version of the TRBA 250 (March 2014 - Chapter 4.2.5 - prevention of needle-stick injuries) promotes the use of safe vascular access devices.

On the following pages, we offer a selection of customizable all-in-one solutions with safety products preventing puncture and cutting injuries for

- central venous catheterization and
- accessory kits.

ProSet provides enhanced user safety through the use of safe instruments.

You can learn more about the directive and the safe instruments that are already part of our product portfolio at our website:

www.safeinfusiontherapy.com



ProSet – we set standards!



ProSet Certofix® Safety

Practice-oriented complete kits for using central venous catheters including safety products have proven to be valuable for a large number of users.




ProSet Accessory Safety

Practice-oriented supplement kits including safety products for central venous catheterization offer a high degree of flexibility.



Central Venous Catheters

Certofix®: Your Benefits



Soft Tip

Made from soft polyurethane, avoids vein perforations and damage to the intima. This means a high level of safety for your patients.

Safsite®

Opens automatically when the Luer connector is attached and closes on disconnection. This reliably prevents the leakage of fluid or blood as well as the entry of air into the system.

Length markings

Printed line-type length markings for determining insertion depth.

Lumina indicators

Printed lumina indicators for accurate identification.

Slide clamps

Safe and easy closing of individual lumina.

Lumen overview

From single lumen catheters for short term intervention to quintuple lumen catheters in the intensive care unit – B. Braun provides solutions for all application areas.



Mono
ø G 14



Duo
ø G 16/16



Trio
ø G 16/18/18



Quattro
ø G 14/18/18/16



Quinto
ø G 16/18/18/18/12

Basic Components of a CVC Kit

Dilator

For expanding tissue structures in the puncture channel.

Safety Seldinger needle


With the proven B. Braun clip technology.

ECG-connection cable

- Determination of catheter placement in real time
- High precision and safety without radiation exposure to patients
- Standard in all Certofix® catheters

Kink-proof guide wire

Extreme flexible and resilient guide wire made from Nitinol with length markings.



Fixation clip

For safe and secure catheter fixation.

Hypodermic Needle-Pro®

Needle with needle protection device.





Safety scalpel

With active safety mechanism.

Guide wire in dispenser

Convenient dispenser for the Seldinger guide wire.

Catheter Sets with Certofix® Mono/Duo/Trio







Contents	Qty
Certofix® – Mono catheter S430 30 cm, ▶ Seldingerkanüle Safety 1.3 x 70 mm, G18 Dilator, guide wire in dispenser, ECG-connection cable, fixation clip, Safety scalpel	1
Extra-large gauze swabs	3
Plum-sized gauze swabs	5
Disposable forceps	1
Fenestrated drape, blue 75 x 90 cm; 9 cm decentral opening; perforated, transparent window 35 x 35 cm	1
Injekt® syringe 10 ml LS	1
Discofix® C-3, blue (included)	1
Silkam® suture, size 3.5=(0) 75 cm, DS30 (curved needle)	1
Examination gown, blue 130 cm	1
Green plastic tray with integrated section	1
Disinfection tray, white	1
Field sheet 75 x 90 cm	1

Contents	Qty
Certofix® – Mono catheter S430 30 cm, ▶ Seldingerkanüle Safety 1.3 x 70 mm, G18 Dilator, guide wire in dispenser, ECG-connection cable, fixation clip, Safety scalpel	1
Gauze compresses 10 x 10 cm	5
Swabsticks	3
Fenestrated drape, blue 100 x 120 cm; 16 x 12 cm decentral opening; perforated, transparent window 35 x 35 cm	1
Omnifix® syringe 10 ml LS	1
Omnifix® syringe 5 ml LL	1
▶ Needle-Pro® 0.7 x 38 mm; G22	1
Sterifix® filter straw 10 cm	1
Discofix® C-3, blue (included)	1
Silkam® suture, size 3.5=(0) 75 cm, DS30 (curved needle)	1
SUSI® needle holder Mayo-Hegar 15 cm	1
Askina® Soft 9 x 10 cm	1
Examination gown, blue 130 cm	1
Green plastic tray with integrated section	1
Disinfection tray, white	1
Field sheet 75 x 90 cm	1

Contents	Qty
Certofix® – Duo catheter S730 30 cm, ▶ Safety Seldinger needle 1.3 x 70 mm, G18 Dilator, guide wire in dispenser, ECG-connection cable, fixation clip, Safety scalpel	1
Gauze compresses 10 x 10 cm	5
Swabsticks	3
Fenestrated drape, blue 100 x 120 cm; 16 x 12 cm decentral opening; perforated, transparent window 35 x 35 cm	1
Omnifix® syringe 10 ml LS	1
Omnifix® syringe 5 ml LL	1
▶ Needle-Pro® 0.7 x 38 mm; G22	1
Sterifix® filter straw 10 cm	1
Discofix® C-, blue (1 lock-fixed, 1 additional)	2
Silkam® suture, size 3.5=(0) 75 cm, DS30 (curved needle)	1
SUSI® needle holder Mayo-Hegar 15 cm	1
Askina® Soft 9 x 10 cm	1
Examination gown, blue 130 cm	1
Green plastic tray with integrated section	1
Disinfection tray, white	1
Field sheet 75 x 90 cm	1

Contents	Qty
Certofix® – Trio catheter S720 20 cm, ▶ Safety Seldinger needle 1.3 x 70 mm, G18 Dilator, guide wire in dispenser, ECG-connection cable, fixation clip, Safety scalpel	1
Gauze compresses 10 x 10 cm	5
Swabsticks	3
Fenestrated drape, blue 75 x 90 cm; 9 cm decentral opening; perforated, transparent window 35 x 35 cm	1
Omnifix® syringe 10 ml LS	1
Omnifix® syringe 5 ml LS	1
▶ Needle-Pro® 0.9 x 38 mm; G20	1
▶ Needle-Pro® 0.7 x 38 mm; G22	1
Discofix® C-3, blue (2 lock-fixed, 1 additional)	3
Dagrofil® suture, size 2/0=(3) 75 cm, DS30 (curved needle)	1
SUSI® needle holder Mayo-Hegar 15 cm	1
Green plastic tray with integrated section	1
Disinfection tray, white	1
Field sheet 75 x 90 cm	1

Catheter Sets with Certofix® Trio






Contents	Qty
Certofix® - Trio catheter S720 20 cm, ▶ Safety Seldinger needle 1.3 x 70 mm, G18 Dilator, guide wire in dispenser, ECG-connection cable, fixation clip, Safety scalpel	1
Gauze compresses 10 x 10 cm	5
Plum-sized gauze swabs	5
Plastic clip	1
2-piece fenestrated drape, blue 75 x 90 cm; central opening 27 x 34 cm	1
Omnifix® syringe 5 ml LS	1
Discofix® C-3, blue (2 lock-fixed, 1 additional)	3
Dagrofil® suture, size 3.5=(0) 75 cm, DS30 (curved needle)	1
SUSI® needle holder Mayo-Hegar 15 cm	1
Askina® Soft 9 x 5 cm	1
Green plastic tray with integrated section	1
Round bowl 60 ml	2
Field sheet 75 x 90 cm	1

Contents	Qty
Certofix® - Trio catheter S720 20 cm, ▶ Safety Seldinger needle 1.3 x 70 mm, G18 Dilator, guide wire in dispenser, ECG-connection cable, fixation clip, Safety scalpel	1
Gauze compresses 10 x 10 cm	5
Swabsticks	3
Fenestrated drape, blue 100 x 120 cm; 16 x 12 cm decentral opening; perforated, transparent window 35 x 35 cm	1
Omnifix® syringe 10 ml LS	1
Omnifix® syringe 5 ml LL	1
▶ Needle-Pro® 0.7 x 38 mm; G22	1
Sterifix® filter straw 10 cm	1
Discofix® C-3, blue (2 lock-fixed, 1 additional)	3
Silkam® suture, size 3.5=(0) 75 cm, DS30 (curved needle)	1
SUSI® needle holder Mayo-Hegar 15 cm	1
Askina® Soft 9 x 10 cm	1
Examination gown, blue 130 cm	1
Green plastic tray with integrated section	1
Disinfection tray, white	1
Field sheet 75 x 90 cm	1

Contents	Qty
Certofix® - Trio catheter S730 30 cm, ▶ Safety Seldinger needle 1.3 x 70 mm, G18 Dilator, guide wire in dispenser, ECG-connection cable, fixation clip, Safety scalpel	1
Gauze compresses 10 x 10 cm	5
Swabsticks	3
Fenestrated drape, blue 100 x 120 cm; 16 x 12 cm decentral opening; perforated, transparent window 35 x 35 cm	1
Omnifix® syringe 10 ml LS	1
Omnifix® syringe 5 ml LL	1
▶ Needle-Pro® 0.7 x 38 mm; G22	1
Sterifix® filter straw 10 cm	1
Discofix® C-3, blue (2 lock-fixed, 1 additional)	3
Silkam® suture, size 3.5=(0) 75 cm, DS30 (curved needle)	1
SUSI® needle holder Mayo-Hegar 15 cm	1
Askina® Soft 9 x 10 cm	1
Examination gown, blue 130 cm	1
Green plastic tray with integrated section	1
Disinfection tray, white	1
Field sheet 75 x 90 cm	1

Contents	Qty
Certofix® - Trio catheter S730 30 cm, ▶ Safety Seldinger needle 1.3 x 70 mm, G18 Dilator, guide wire in dispenser, ECG-connection cable, fixation clip, Safety scalpel	1
Extra-large gauze swabs	3
Plum-sized gauze swabs	5
Disposable forceps	1
Fenestrated drape, blue 75 x 90 cm; 9 cm decentral opening; perforated, transparent window 35 x 35 cm	1
Injekt® syringe 10 ml LS	1
Discofix® C-3, blue (included)	3
Safsite® safety connector (lock-fixed)	2
Silkam® suture, size 3.5=(0) 75 cm, DS30 (curved needle)	1
Examination gown, blue 130 cm	1
Green plastic tray with integrated section	1
Disinfection tray, white	1
Field sheet 75 x 90 cm	1

Catheter Sets with Certofix® Quattro/Quinto




Contents	Qty
Certofix® – Quattro catheter S820 20 cm, ▶ Safety Seldinger needle 1.3 x 70 mm, G18 Dilator, guide wire in dispenser, ECG-connection cable, fixation clip, Safety scalpel	1
Extra-large gauze swabs	3
Plum-sized gauze swabs	5
Disposable forceps	1
Fenestrated drape, blue 75 x 90 cm; 9 cm decentral opening; perforated, transparent window 35 x 35 cm	1
Injekt® syringe 10 ml LS	1
Discofix® C-3, blue (included)	4
Safsite® safety connector (lock-fixed)	3
Silkam® suture, size 3.5=(0) 75 cm, DS30 (curved needle)	1
Examination gown, blue 130 cm	1
Green plastic tray with integrated section	1
Disinfection tray, white	1
Field sheet 75 x 90 cm	1

Contents	Qty
Certofix® – Quinto catheter S1220 20 cm, ▶ Safety Seldinger needle 1.3 x 70 mm, G18 Dilator, guide wire in dispenser, ECG-connection cable, fixation clip, Safety scalpel	1
Gauze compresses 10 x 10 cm	5
Plum-sized gauze swabs	5
Plastic clamp	1
2-piece fenestrated drape, blue 75 x 90 cm; central opening 27 x 34 cm	1
Omnifix® syringe 5 ml LS	1
Discofix® C-3, blue (4 lock-fixed, 1 additional)	5
Dagrofil® suture, size 3.5=(0) 75 cm, DS30 (curved needle)	1
SUSI® needle holder Mayo-Hegar 15 cm	1
Askina® Soft 9 x 5 cm	1
Green plastic tray with integrated section	1
Round bowl 60 ml	2
Field sheet 75 x 90 cm	1


Contents	Qty
Certofix® – Quinto catheter S1220 20 cm, ▶ Safety Seldinger needle 1.3 x 70 mm, G18 Dilator, guide wire in dispenser, ECG-connection cable, fixation clip, Safety scalpel	1
Gauze compresses 10 x 10 cm	5
Swabsticks	3
Fenestrated drape, blue 75 x 90 cm; 9 cm decentral opening; perforated, transparent window 35 x 35 cm	1
Omnifix® syringe 10 ml LS	1
Omnifix® syringe 5 ml LS	1
▶ Needle-Pro® 0.9 x 38 mm; G20	1
▶ Needle-Pro® 0.7 x 38 mm; G22	1
Discofix® C-3, blue (4 lock-fixed, 1 additional)	5
Dagrofil® suture, size 2/3=(3) 75 cm, DS30 (curved needle)	1
SUSI® needle holder Mayo-Hegar 15 cm	1
Green plastic tray with integrated section	1
Disinfection tray, white	1
Field sheet 75 x 90 cm	1

ProSet Accessory Safety

Accessory Kits for Central Venous Catheters



Article number: 4897724
Sales unit: 10 items




Article number: 4897725S
Sales unit: 10 items


Contents	Qty
▶ Safety scalpel	1
Gauze compresses 10 x 10 cm	5
Plum-sized gauze swabs	5
Plastic clamp	1
2-piece fenestrated drape, blue 75 x 90 cm; central opening 27 x 34 cm	1
Omnifix® syringe 5 ml LS	1
Dagrofil® suture, size 3.5=(0) 75 cm, DS30 (curved needle)	1
SUSI® needle holder Mayo-Hegar 15 cm	1
Askina® Soft 9 x 5 cm	1
Green plastic tray with integrated section	1
Round bowl 60 ml	2
Field sheet 75 x 90 cm	1

Contents	Qty
▶ Safety scalpel	1
Gauze compresses 10 x 10 cm	5
Plum-sized gauze swabs	5
Plastic clamp	1
Fenestrated drape, blue 100 x 120 cm; 16 x 12 cm decentral opening; perforated, transparent window 35 x 35 cm	1
Omnifix® syringe 5 ml LS	1
Dagrofil® suture, size 3.5=(0) 75 cm, DS30 (curved needle)	1
SUSI® needle holder Mayo-Hegar 15 cm	1
Askina® Soft 9 x 5 cm	1
Green plastic tray with integrated section	1
Round bowl 60 ml	2
Field sheet 75 x 90 cm	1


Accessory Kits for Central Venous Catheters




Article number: 4890002S
Sales unit: 10 items



Article number: 4895696S
Sales unit: 8 items



Article number: 4892000
Sales unit: 10 items



Article number: 4892001
Sales unit: 10 items

Contents	Qty
Gauze compresses 10 x 10 cm	5
Non-woven slit compresses 10 x 10 cm	5
Plum-sized gauze swabs	5
Plastic clamp	1
Fenestrated drape, blue 75 x 90 cm; 9 cm central opening; transparent window 35 x 35 cm	1
Omnifix® syringe 10 ml LS	3
Omnifix® syringe 5 ml LS	1
Omnifix® syringe 5 ml LL	1
▶ Needle-Pro® 0.9 x 38 mm; G20	1
▶ Needle-Pro® 0.7 x 38 mm; G22	1
Discofix® C-3, blue	3
Dagrofil® suture, size 3=2/0 75 cm, DS24 (curved needle)	1
Askina® Soft 9 x 5 cm	1
Examination gown, blue 130 cm	1
Green plastic tray with integrated section	1
Disinfection tray, white	1
Field sheet 75 x 90 cm	1

Contents	Qty
▶ Safety scalpel	1
Gauze compresses 10 x 10 cm	10
Swabsticks	3
Fenestrated drape, blue 75 x 90 cm; 9 cm central opening; transparent window 35 x 35 cm	1
Injekt® syringe 10 ml LS	1
Injekt® syringe 5 ml LS	1
Injekt® syringe 2 ml LS	1
▶ Needle-Pro® 0.5 x 38 mm; G25	1
Discofix® C-3, blue	2
Discofix® C-3, red	1
Needle-Block-Foam	1
Askina® Soft 9 x 10 cm	1
Examination gown, blue 150 cm	1
Green plastic tray with integrated section	1
Disinfection tray, white	1
Field sheet 75 x 90 cm	1

Contents	Qty
Gauze compresses 7.5 x 7.5 cm	5
Swabsticks	3
Fenestrated drape, blue 70 x 90 cm; 7 x 14 cm oval decentral opening	1
Omnifix® syringe 20 ml LS	1
Omnifix® syringe 10 ml LS	2
Omnifix® syringe 5 ml LS	1
▶ Needle-Pro® 0.8 x 38 mm; G21	1
▶ Needle-Pro® 0.5 x 25 mm; G25	2
Sterifix® filter straw 10 cm	1
Ultrasound cover 13 cm x 120 cm	1
Green plastic tray with integrated section	1
Disinfection tray, white	1
Field sheet 75 x 90 cm	1

Contents	Qty
Gauze compresses 7.5 x 7.5 cm	5
Swabsticks	3
Fenestrated drape, blue 70 x 90 cm; 7 x 14 cm oval decentral opening	1
Omnifix® syringe 20 ml LS	1
Omnifix® syringe 10 ml LS	2
Omnifix® syringe 5 ml LS	1
▶ Needle-Pro® 0.8 x 38 mm; G21	1
▶ Needle-Pro® 0.5 x 25 mm; G25	2
Sterifix® filter straw 10 cm	1
Ultrasound cover 13 cm x 120 cm, Ultrasound gel 20 g	1
Askina® Derm 10 x 12 cm	1
Green plastic tray with integrated section	1
Disinfection tray, white	1
Field sheet 75 x 90 cm	1

B. Braun Melsungen AG | Hospital Care | 34209 Melsungen | Germany
Tel. +49 5661 71-0 | www.bbraun.com



Learn more about ProSet online – Scan the QR code
with your mobile device or visit our website:
www.proset-bbraun.com