

# Askina® Soft | Askina® Soft Clear

## Sterile post-operative wound dressing



Wound Management

# Askina® Soft | Askina® Soft Clear

Sterile post-operative wound dressing



**Askina® Soft** is a sterile post-operative dressing combining a breathable fleece material and a non-adherent absorptive pad.

### Indications

Management of post-operative wounds, and all sorts of acute wounds such as cuts, lacerations, sutures, abrasions, ...


- Elastic wound dressing
- Breathable fleece material
- Skin friendly adhesive
- Non-adherent wound contact layer
- Easy to apply, even on difficult locations


**Askina® Soft Clear** is a sterile post-operative dressing combining a semi-permeable transparent film and a non-adherent absorptive pad.

### Indications

Management of post-operative wounds, and all sorts of acute wounds requiring a waterproof barrier and a visual control.

- Easy to apply, even on difficult locations
- Skin friendly acrylic adhesive
- Allows the patient to shower with the dressing in situ
- Bacteria proof barrier
- Allows visual control

Askina® Soft	Size	Pcs/Pack	Article no
	7.5 cm x 5 cm	50	9086480
	9 cm x 5 cm	50	9086501
	9 cm x 10 cm	50	9086510
	9 cm x 15 cm	40	9086528
	9 cm x 20 cm	30	9086536
	9 cm x 25 cm	30	9086544
	9 cm x 30 cm	30	9086552

Askina® Soft Clear	Size	Pcs/Pack	Article no
	7.5 cm x 5 cm	50	9086803
	9 cm x 5 cm	50	9086812
	9 cm x 10 cm	50	9086821
	9 cm x 15 cm	40	9086830
	9 cm x 20 cm	30	9086839
	9 cm x 25 cm	30	9086848
	9 cm x 30 cm	30	9086857