

# Askina® Soft I.V. | Askina® Soft Clear I.V. I.V. cannula fixation dressing



Wound Management

# Askina® Soft I.V. | Askina® Soft Clear I.V.

## I.V. cannula fixation dressing

Skin friendly adhesive

Slit to wrap around tubing

Extra pad for cushioning

Non-adherent absorptive pad  
 – The puncture site stays dry  
 – Prevents puncture site adherence

Ideal for fixation of Vasofix® Safety



### Askina® Soft I.V.




Sterile I.V. cannula fixation dressing, with an integrated non-adherent absorptive pad. Made of soft conformable non-woven material.


### Askina® Soft Clear I.V.



Sterile transparent version of Askina® Soft I.V., for an alternative use.

- Bacteria proof barrier
- Visual inspection of the puncture area
- Breathable material

Askina® Soft I.V.	Size	Pcs/Pack	Article no
	8 cm x 6 cm	50	9086919

Askina® Soft Clear I.V.	Size	Pcs/Pack	Article no
	8 cm x 6 cm	50	9086893