



VascuGraft® FLOW

Heparin coated vascular prostheses. Your right choice for all relevant indications.

Indications

- Peripheral arterial reconstructions
- Extra-anatomical procedures
- Dialysis shunt placement

VascuGraft® FLOW

Increased graft patency due to a thromboresistant inner luminal surface

VascuGraft FLOW is characterized by an innovative heparin coating on the luminal side of the prostheses: The covalently bound heparin on the inner surface of VascuGraft FLOW enhances the absorption of specific plasma proteins and thus enables the fast formation of a natural vessel-like anti-thrombogenic layer.

As further important function of heparin, the formation of the procoagulant thrombin as contributing factor to thrombus formation will be inhibited.

