





SOLUTION FOR PATIENTS WITH HIGH OUTPUT Quality of Life our main asset

SOLUTION FOR PATIENTS WITH HIGH OUTPUT

Your patients are looking for:

- Good sleep at night
- Better recovery and QoL
- Security and freedom

Did you know that patients with a High Output Stoma have to empty their pouch 8 to 25 times per 24h?*

Made by nurses – Approved by patients

85.7% of patients feel more secure.**

21/24 patients use it at night & for 93.7% of patients enhanced their Quality of Life. **



66.6% enhanced their autonomy & for 90% of patients it provides optimum freedom of movements! **

* Slater R High-output stomas : challenges with a large laparostomy wound. BJN 2011; 21(6):S26–S33 Calculation based on an output from 2L to 8L with an emptying every 300 to 600ml. ** Report on Post Market User Evaluation Flow Collector, July 2014 – Data on file

Making high output patients management easier!

- + Secure connection between Flexima® pouch and Flow Collector
- + Enhanced autonomy thanks to the 2L distal collecting pouch

Flow Collector

+ Discrete

- Adjustable (every 30cm) white opaque tube
- White opaque backing

+ Easy monitoring

• Clear printed graduate scale

• Indications every 250mL - Up to 2L approx. (72h collection)

+ Practical & hygienic outlet • Outlet located at the bottom end:

- - 2 grips to facilitate the opening
 - System to hold the cap when emptying
- Individually wrapped for an easy handling in outpatient line

Specific outlet w	ith notches
to ensure strong	connection

FLOW COLLECTOR

100% of patients are

satisfied with collecting

capacity**

Up to 1.5m tube

arphi $\dot{\sim}$

FLOW COLLECTO

BBRAUN

85% of patients are satisfied with the emptying system *

BIBB

Compatible with High Flow stoma pouches :



Both versions are compatible with the Flow Collector

+ Secure & hygienic outlet:

- Atraumatic edges thanks to the soft material
- Compact shape of the outlet for better comfort
- Secure connection to the Flow Collector

Large (50mm) & direct flow for an easy emptying

Large 13.3 mm diameter to evacuate debris -

SOLUTION FOR PATIENTS WITH HIGH OUTPUT

			NEW	
High Flow pouches (Box of 30)	Diameter	Transparent	Split Cover Beige w/ inspection window	Beige
Flexima [®] High Flow	15-75	^O 42971	42970	-
Flexima [®] 3S High Flow	55	^O 933155	933055	-
Flexima [®] 3S High Flow	65	^O 933165	933065	-
Flexima [®] 3S High Flow	80	^O 933180	933080	-
Flexima® Key Maxi Drainable	50	^O 62650	-	62550
Flexima® Key Maxi Drainable	60	^O 62660	-	62560
Flexima® Key Maxi Drainable	80	^O 62680	-	-

Flexima® 3S Base Plates*	Diameter	References	Flexima® Key Base Plates*	Diameter	References
Flexima® 3S Flat	45/ 15-30	936415	Flexima® Key flat	40/ 15-35	62040
Flexima® 3S flat	55/ 15-40	936515	Flexima® Key flat	50/ 15-45	62050
Flexima® 3S flat	65/ 15-50	936615	Flexima® Key flat	60/ 15-55	62060
Flexima® 3S flat	80/ 15-65	936815	Flexima® Key flat	80/ 15-75	62080
Flexima® 3S convex	45/ 15-26	937415	Flexima [®] Key convex	50/ 15-25	62052
Flexima [®] 3S convex	55/ 15-36	937515	Flexima [®] Key convex	60/ 15-35	62062
Flexima® 3S convex	65/ 15-46	937615	Flexima [®] Key convex	60/ 15-45	62063
Flat: box of 10/ Convex: box of 5			*Pre-cut base plates are also ava	ilable / Flat & Conve	x: box of 5

Flow Collector	Packaging	References
	Box of 5	039901
	Box of 15	039900

This document, its contents, including institutional data, information, trademarks and logos mentioned herein are the exclusive property of B. Braun. Any representation and / or reproduction, total or partial, of this document and its contents without the express prior consent of B. Braun, is strictly prohibited and constitutes an infringement of the intellectual property rights of B. Braun. Non binding documents and photographs. For healthcare professional use only. Please read carefully the instructions for use.

Manufactured and distributed by:

B. Braun Medical | 204 avenue du Maréchal Juin | 92100 Boulogne-Billancourt | FRANCE Tél. 01 41 10 53 00 | Fax 01 41 10 53 99 | www.bbraun.com B. Braun Medical | SAS with a share capital of 31 000 000 € | RCS Nanterre 562050856

For most countries a specific letter is added at the end of the radical reference. Availability may vary from one country to another.

All volumes are measured acording to ISO 8670-2 standard. This medical device is a regulated health care product and on the basis of this regulation is certified CE.

The products commercialized by B. BRAUN MEDICAL shall exclusively be used in accordance with the package inserts which have been updated and included in the boxes.